

Immuno-Assay Validation Report

Facility:
N°

Study director:

BIOMARKER / ANALYTE: *Precise the most common names of the analyte which will be measured by the assay.*

TYPE OF ASSAY:

- Microtiter-competitive ELISA
 Microtiter-sandwich ELISA
 Bead-based immunoassay
 Microtiter multiplex MSD

TYPE OF VALIDATION:

- Pre-exploratory validation *(for PoT studies)*
 Exploratory validation *(during PoT studies)*
 Pre-confirmatory validation *(for PoP studies)*

REAGENTS:

- Kit.....Source:.....Batch number:.....Date of production:.....
 Antibody.....Source:.....Batch number:.....Date of production:.....
 Antigene.....Source:.....Batch number:.....Date of production:.....
 Standard.....Source:.....Batch number:.....Date of production:.....

Date of validation study completed:

Operator(s):

Date:

Validation procedure

See SOP: SAFE-T Standard Validation Procedure for Biomarker immunoassays.

Objectives / Intended use of the assay:

*To demonstrate that the performance characteristics of the assay make it fit for the intended analytical application.
 Define the intended use of the immune-assay (biological fluid, disease, etc).*

Immuno-assay validation plan

Related to Immuno-Assay Validation Plan N°:

Immuno-Assay validation plan written by:

Validation results

Attach all data record forms and charts (table of all analytical runs with analysis dates).

Validation standards:

Description of the validation standards and matrix used for validation study.

The source and lot number, expiration date, certificates of analyses when available, and/or internally or externally generated evidence of identity and purity should be furnished for each reference standard.

Calibration curve

The calibration curve should be shown with the fitting model clearly indicated (e.g. 4PL, 5PL, etc...).

A table for observed and fitted data should be provided.

Assay dynamic range

LLOQ:

ULOQ:

*Give a short reasoning for the determination of the LLOQ and ULOQ.***Sensitivity:**

Limit of detection (LOD):

*(minimum detectable dose).**Give a short overview of procedure used to determine the LOD.***Precision****Intra-assay precision:***The CV(%) of the middle concentration should be given.***Inter-assay precision:***A table with the mean value, the SD and CV for 5 samples with 3 different concentrations should be shown for inter-assay precision.***Operator-to-operator variability:***Show the effect of the operator on inter-assay variation with at least 2 operators***Recovery: %***Mean recovery from triplicate measurements for 3 different concentrations covering the assay range***Selectivity / Specificity: % error***Potential interfering substance should be cited.***Parallelism***The parallelism curves and results should be shown.***Dilutional linearity***The dilutional linearity curve and results should be shown.***Robustness****- Stability**

- *Freeze and Thaw Stability:*
- *Short-Term Temperature Stability:*
- *Long-Term Stability:*
- *Stock Solution Stability:*

- Reliability*If data are available***Suitable matrix for analyte detection***Give a short reasoning for suitable matrix types.***Validation evaluation: Acceptance criteria pass/fail***Compare results to acceptance criteria.*

Criteria	Results	Pass/Fail
(See validation plan)		
Calibration curve: appropriate model _____		<input type="checkbox"/> / <input type="checkbox"/>
Assay dynamic range: Fitted to target range _____		<input type="checkbox"/> / <input type="checkbox"/>
LLOQ: lowest calibrator concentration of the back calculated concentrations that does not exceed a CV of 20% while the recovery is within 75-125% _____		<input type="checkbox"/> / <input type="checkbox"/>
ULOQ: highest calibrator concentration of the back calculated concentrations that does not exceed a CV of 20% while the recovery is within 75-125% _____		<input type="checkbox"/> / <input type="checkbox"/>
Sensitivity: LOD = 3 times the standard deviation of the background _____		<input type="checkbox"/> / <input type="checkbox"/>
Accuracy: _____		<input type="checkbox"/> / <input type="checkbox"/>
Precision:		
Intra-assay precision: 15% CV for a sample in midrange _____		<input type="checkbox"/> / <input type="checkbox"/>
Inter-assay precision: 20% CV for the VS samples except for the LLOQ < 30% _____		<input type="checkbox"/> / <input type="checkbox"/>
Recovery: 100±20% _____		<input type="checkbox"/> / <input type="checkbox"/>
Selectivity / Specificity: %error within reasonable range _____		<input type="checkbox"/> / <input type="checkbox"/>
Parallelism: _____		<input type="checkbox"/> / <input type="checkbox"/>
Dilution linearity: _____		
Stability:		
Freeze and Thaw Stability: deviation should be within ±25% _____		<input type="checkbox"/> / <input type="checkbox"/>
Short-Term Temperature Stability: deviation should be within 15% _____		<input type="checkbox"/> / <input type="checkbox"/>
Long-Term Stability: _____		

deviation should be within $\pm 25\%$ _____ <input type="checkbox"/> / <input type="checkbox"/>
Stock Solution Stability: deviation should be within $\pm 15\%$ _____ <input type="checkbox"/> / <input type="checkbox"/>
Validation conclusion <i>Is the assay validated?</i>
Protocol deviations / descriptions and justifications
Additional data and information relevant to the validation study <i>All additional data needed for the validation report (Information about antibodies or reagent, Assay procedure, Calculations and Statistical Analyses, ...)</i>
Validation report approval Study director Name: Facility: Validation report approved: <input type="checkbox"/> Yes <input type="checkbox"/> No Comments: Date: Study director signature: